











## Imported Baxter 1000mL 0.9% Sodium Chloride in Viaflex (MPC: 2B1324X)

Replacing: Sodium Chloride 0.9% w/v intravenous infusion BP (MPC: FKE1324)

Overpouch opening			
			
Using aseptic technique: Hold the overpouched bag vertically. Grasp the top seal of the overpouch close to one of the slits	Tear the overpouch from the slit downwards in a smooth movement		Shake the exposed bag from the overpouch



Spike connection				
				
Using aseptic technique, remove the protective cap from the IV administration set and close the roller clamp	Remove the blue protective cap, by pulling it downwards until it pops off	Hold the administration port on the collar so as not to occlude the administration port tubing	Firmly insert the spike of the set straight into the administration port. You may hear or feel a click as the inner membrane is punctured.  <b>Important: Do not twist the set during insertion, but do twist the set to aid set removal from bag.</b>	

Baxter Healthcare Ltd  
Wallingford Road  
Compton  
Newbury  
Berkshire  
RG20 7QW

Adverse Events and any suspected defective medicines should be reported. Reporting forms and information can be found at [www.mhra.gov.uk/yellowcard](http://www.mhra.gov.uk/yellowcard). Adverse Events relating to Baxter products can also be reported direct to Baxter Pharmacovigilance on +44 (0)1635 206360, or by email to [vigilanceuk@baxter.com](mailto:vigilanceuk@baxter.com). Any drug product quality complaints (including suspected defective medicines) relating to Baxter products can be reported directly to the Baxter Country Quality Assurance Team on +44 (0)1604 704603, or by email to [UK\\_SHS\\_QA\\_Complaints@baxter.com](mailto:UK_SHS_QA_Complaints@baxter.com). Alternatively, please report directly to your Baxter Representative, who will take the details and forward to the Baxter Country Quality Assurance Team.

To download Product Leaflet and Comparison Guide please refer to the Baxter Product Catalogue: <https://ecatalog.baxter.com/ecatalog/home>. Search for Product Code FKE1324; the information can be found within the Product Characteristics section, under Clinical Data and FAQ.